

New



Yes, you can.

➤ *Invacare*<sup>®</sup> **Action**<sup>®</sup> **2NG**

The Next Generation has arrived!

The all-new **Action 2NG** has been developed from the tried-and-tested Action 3NG concept. The specifications meet the basic needs of the users whilst refusing to compromise on safety or comfort. The next generation **Action 2NG** shares the same essential platform as the Action 3NG, offering a low total weight for easy handling combined with good driving abilities. Due to the shared use of common components, the **Action 2NG** perfectly fits into the Action family.



## Light and comfortable

With a weight of just 15.kg, the Invacare® **Action 2NG** fulfills all essential functional requirements. The stable folding frame means the wheelchair is easy to handle, and the lite weight makes transportation trouble-free. Cross-compatibility with other members of the Action family ensures that the **Action 2NG** can be easily adapted to meet individual user needs.

## Quality of life

The tried-and-tested technology of the Action family products allow the user to go wherever they desire.

## Features



**Swing in and out footrest**  
The 80° footrest swings in and out for easy transfers in or out of the chair.



**Adjustable center of gravity**  
The newly-designed rear wheel plate can be adjusted by 5 cm.



**Anti-tippers**  
Provide extra safety on slopes or uneven ground.



## Technical Data

15", 16", 17", 18", 19"	16 - 17"	48.5, 51 cm	43 cm
20 - 26 cm	35 - 48 cm	80° / 0° - 80°	7° / 0 - 30°
SW+ 20.5 cm	93.5 cm	82 / 84.5 cm	41 cm
15kg	max. 125 kg	Crash tested ISO 7176-19*	

\* Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.

## Frame Colour



Anthracite Black



Pearl Grey

Please note colors may vary slightly from those displayed above.

## Upholstery



Black Nylon

## Your Local Dealer:



Yes, you can.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Invacare Australia Pty Ltd ABN 45 074 676 378.  
PO Box 5002, 1 Lenton Place, North Rocks, NSW 2151, Australia.  
Freephone: 1800 460 460, Fax: 02 8839 5353,  
Email: sales@invacare.com.au, www.invacare.com.au

Invacare New Zealand,  
PO Box 62-124, 4 Westfield Place, Mt Wellington, Auckland, New Zealand.  
Freephone: 0508 468 222, Freefax: 0508 807 788  
Email: sales@invacare.co.nz, www.invacare.co.nz

Invacare Asia Ltd,  
PO Box 5002, 1 Lenton Place, North Rocks, NSW 2151, Australia.  
Phone: +61 2 8839 5330, Fax +61 2 8839 5331,  
Email: asiasales@invacare.com, www.invacare.com